

GP DETERGENT

GREEN NEUTRAL WASHING-UP LIQUID

CLOVER 
CHEMICALS LTD
'Effective Cleaning Chemicals'
www.cloverchemicals.com

- Mild - ideal for hand use
- Contains additives for effective grease removal
- Long lasting foam

USES

A detergent concentrate containing additives for efficient grease removal.


A thick, green liquid consisting of detergents, grease-dissolving surfactants, foam boosters, preservative and dye. GP detergent is pH-balanced for mildness to hands.

Designed for cleaning glassware, cutlery, crockery, utensils and worktops, GP Detergent is ideal for use on aluminium, glass, stainless steel, vinyl, plastic, ceramic, wood and enamel.

DIRECTIONS

DISHWASHING Dilute 1 part GP Detergent with between 750 to 1000 parts warm water in a washing-up bowl. Add items to be washed and allow a few minutes to soak. Clean items using a sponge, brush or pad.



GENERAL CLEANING - Dilute 1 part GP Detergent with 300 parts cold or warm water.



DILUTIONS

GP DETERGENT : WATER	
Dishwashing 1:750 - 1:1000	General Cleaning 1:300

Code: 418



PACK SIZES

